

HOOC-PEG-COOH (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HOOC-PEG-COOH (MW 40000) is a compound with two active carboxyl groups capable of selectively reacting with amino groups in peptide coupling agents. It serves as an excellent crosslinker for the PEGylation of proteins and peptides, as well as for the modification of nanoparticles and surfaces.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481