

DSPE-PEG3400-THRPPMWSPVWP

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG3400-THRPPMWSPVWP is a PEG compound made up of DSPE and the transferrin receptor (TfR) peptide. The sequence THRPPMWSPVWP binds to TfR and is then internalized into cells that express TfR.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481