

DSPE-PEG3400-BR2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3400-BR2 is a PEG compound comprised of DSPE and BR2. BR2 is a targeting peptide made up of 17 amino acids with the sequence RAGLQFPVGRLRRLLR. It possesses the capability to assist cell membrane penetration and can specifically target certain cells or tissues.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481