

DSPE-PEG3400-iRGD

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3400-iRGD is a PEG compound composed of DSPE and the α v-integrin targeting peptide (iRGD). The iRGD peptide binds to α v-integrin, undergoing proteolytic cleavage within tumors to produce CRGDK/R, which interacts with neuropilin-1, thereby exhibiting tumor-targeting and penetrating properties. DSPE-PEG3400-iRGD is applicable for drug delivery.
Targets(IC50)	Integrin,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481