

DSPE-PEG-NHS ammonium (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-NHS (ammonium) (MW 10000) is a widely used reactive phospholipid commonly employed for conjugating antibodies, proteins, peptides, or other substrates to the surface of liposomes and other lipid-PEG nanoparticles.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481