

SjDX5-53 control peptide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SjDX5-53 control peptide serves as a control for the TLR2 activator SjDX5-53, without targeting specific receptors or influencing regulatory T cells (Tregs) or dendritic cells. It is useful for excluding nonspecific effects and confirming SjDX5-53's specific role in Tregs induction and autoimmune disease mitigation. This peptide is applied in research on autoimmune diseases like colitis and psoriasis.
Targets(IC50)	TLR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481