

Mal-PEG-NH2 TFA (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mal-PEG-NH2 TFA (MW 40000) is a linear heterobifunctional PEG compound featuring maleimide and ?NH2 groups. The maleimide group can react with thiol, SH, mercapto, or sulfhydryl groups, while the amine group is reactive with carboxylic acids or activated NHS esters.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481