

4-Arm PEG-SH (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Arm PEG-SH (MW 40000) (4-Arm PEG-thiol (MW 40000)) is a multi-arm PEG derivative featuring thiol groups at the end of its four arms. The reactive free thiol, SH or sulfhydryl groups, selectively interact with maleimide and surfaces of transition metals such as gold and silver.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481