

Alkyne-PEG-NH2 (MW 20000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alkyne-PEG-NH2 (MW 20000) is an alkyne-functionalized, bifunctional PEG derivative used for the modification of proteins, peptides, and other materials. The alkyne group can react with azides in a copper-catalyzed aqueous reaction. Pegylation can enhance solubility and stability while reducing the immunogenicity of peptides and proteins.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481