

## Silane-PEG-COOH (MW 2000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Silane-PEG-COOH (MW 2000) is applicable for modifying proteins, peptides, and other materials with active groups. The carboxyl group (-COOH) readily forms stable amide bonds with amino groups and ester bonds with hydroxyl groups. Silane-functionalized PEG (Silane-PEG-X) is also suitable for modifying surfaces like glass and silicon.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481