

mPEG-CM (MW 40000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	mPEG-CM (MW 40000) [mPEG-COOH (MW 40000)] terminates with a reactive carboxyl group (-COOH) that can form stable amide bonds with amino groups or ester bonds with hydroxyl groups. mPEG-CM (MW 40000) is suitable for drug delivery applications.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481