

DSPE-PEG3400-T7

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG3400-T7 is a PEGylated compound composed of DSPE and the transferrin receptor peptide (T7). The T7 peptide (HAIYPRH) specifically binds to the transferrin receptor (TfR). DSPE-PEG3400-T7 is utilized for drug delivery applications.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481