

DSPE-PEG3400-LTLRWVGLMS

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | DSPE-PEG3400-LTLRWVGLMS is a PEG compound composed of DSPE and the decapeptide (LTLRWVGLMS). The receptor for LLRWVGLMS is chondroitin sulfate proteoglycan NG2. LLRWVGLMS demonstrates homing to pericytes associated with tumor vasculature. DSPE-PEG3400-LTLRWVGLMS can be utilized for drug delivery. |
| Targets(IC50) | Liposome  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481