

4-Arm-PEG-FA (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Arm-PEG-FA (MW 3400) (4-Arm-PEG-folate (MW 3400)) is a PEG derivative where each PEG chain terminates with a Folic acid molecule. Folic acid has a high affinity for folate receptors, which are utilized for targeted drug delivery across cell membrane receptors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481