

HOOC-PEG-SH (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HOOC-PEG-SH (MW 10000) [HOOC-PEG-Thiol (MW 10000)] is a functional thiol PEG derivative featuring terminal carboxyl groups. The carboxyl group can react with amines or hydroxyls to form stable amide bonds or less stable ester bonds. The PEG linkage between the thiol and carboxyl groups offers excellent water solubility, flexible linker distance, and enhanced stability.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481