

HOOC-PEG-COOH (MW 3400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HOOC-PEG-COOH (MW 3400) possesses two reactive carboxyl groups at its ends, enabling selective reaction with amino groups in peptide coupling agents. This compound serves as an excellent crosslinker, suitable for the PEGylation of proteins and peptides, as well as for nanoparticle and surface modification.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481