

Boc-NH-PEG-COOH (MW 2000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | Boc-NH-PEG-COOH (MW 2000) is a linear heterobifunctional PEG compound featuring a Boc-protected amine and a carboxylic acid group. This compound serves as a significant crosslinking agent featuring a PEG chain. |
| Targets(IC50) | Others   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481