

Alkyne-PEG-COOH (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alkyne-PEG-COOH (MW 10000) serves as a heterobifunctional PEG crosslinker, designed for attaching functional PEG to biomolecules, particles, and other material surfaces. The carboxyl group effectively reacts with amine groups to create stable amide bonds.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481