

mPEG-amine (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mPEG-amine (MW 40000) (mPEG-NH ₂ (MW 40000)) is a monofunctional polyethylene glycol with a terminal amine group. It forms stable amide bonds when coupled with NHS esters. Additionally, in the presence of an activator, the amine group can react with carboxyl groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481