

Fmoc-NH-PEG-COOH (MW 2000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Fmoc-NH-PEG-COOH (MW 2000) is a linear heterobifunctional PEG compound featuring an amine protected at both ends by Fmoc and a carboxylic acid. It serves as a significant crosslinker with a PEG chain. The protected amine can be regenerated under mild basic conditions.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481