

Biotin-PEG-COOH (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-PEG-COOH (MW 20000) is a linear heterobifunctional PEG compound featuring biotin and carboxylic acid groups. Biotin binds with high specificity and affinity to avidin and streptavidin. In the presence of a coupling agent, Biotin-PEG-COOH (MW 20000) can biotinylate antibodies, proteins, and other macromolecules containing primary amines.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481