

Alkyne-PEG-COOH (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Alkyne-PEG-COOH (MW 2000) is a heterobifunctional PEG crosslinker that enables the attachment of functional PEG chains to biomolecules, particles, and various material surfaces. The carboxyl group can react with amino groups to form stable amide bonds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481