

mPEG-DMPE (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mPEG-DMPE (MW 1000) is a polyethylene glycol-conjugated lipid featuring a methyl group at the opposite end of the PEG chain. It is utilized in the creation of lipid nanoparticles or liposomes for targeted drug delivery applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481