

## 6-Arm-PEG-Mal (MW 1000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

|               |   |
|---------------|---|
| Description   | 6-Arm-PEG-Mal (MW 1000) is a PEG derivative modified with Maleimide, which forms stable thioether bonds with sulfhydryl groups (-SH). PEG modification enhances solubility and stability, while reducing the immunogenicity of proteins and peptides. |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481