

Chitosan ( $\geq 98\%$  deacetylated, mushroom)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder:  $-20^{\circ}\text{C}$  for 3 years | In solvent:  $-80^{\circ}\text{C}$  for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Chitosan ( $\geq 98\%$ deacetylated, mushroom) is a biochemical reagent that can be used as a flocculant, protein precipitant, encapsulating agent, and aqueous thickener.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481