

Endo-BCN-PEG4-Pomalidomide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Endo-BCN-PEG4-Pomalidomide is a PEG derivative comprised of BCN, PEG4, and Pomalidomide. The BCN group exhibits high reactivity with azide-labeled molecules. Its amino group reacts with carboxylic acids, activated NHS esters, and carbonyl compounds (ketones, aldehydes), among others.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481