

T4 Endonuclease V

Chemical Properties

CAS No. : 52227-85-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	T4EndonucleaseV is a bifunctional DNA glycosylase that exhibits both DNA N-glycosylase and AP lyase activities. It is capable of repairing cyclobutane pyrimidine dimers induced by ultraviolet (UV) radiation in DNA.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481