

CP-EPS8-NLS

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CP-EPS8-NLS is a peptide with anti-AML properties. It disrupts ESP8-related signaling and exhibits anti-AML activity across various AML cell types. Additionally, CP-EPS8-NLS demonstrates potent antitumor effects in xenograft tumor models. It effectively downregulates EPS8 expression along with its downstream pathways.
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481