

4-Arm-PEG-Silane (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Arm-PEG-Silane (MW 20000) is a silane-modified PEG derivative used for surface modification of glass, silicon, or other materials through the reaction of hydroxyl groups with ethoxy (methoxy) silane. PEG effectively inhibits nonspecific surface binding.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481