

## C5 Protein (human)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	C5 Protein (human) is a naturally glycosylated polypeptide composed of two disulfide-linked chains. It is essential for the formation of the membrane attack complex (MAC) and is activated via all three pathways of complement activation.
Targets(IC50)	Complement System

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481