

6-Arm-PEG-Silane (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	6-Arm-PEG-Silane (MW 10000) is a silane-modified PEG derivative that can react with hydroxyl groups and ethoxy (methoxy) silanes to modify surfaces such as glass, silicon, and others. PEG functions to inhibit non-specific binding on surfaces.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481