

## Lenbelintide

## Chemical Properties

CAS No. : 2793341-99-6

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | Lenbelintide is a dual amylin/calcitonin receptor agonist, useful in studies related to obesity. |
| Targets(IC50) | Amylin Receptor  |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481