

SA-21

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SA-21 is a polymeric peptide consisting of 21 amino acids (Ac-RLIEDICLPRWGCLWEDD-NH <sub>2</sub> ) capable of binding to the ligand of urokinase-type plasminogen activator, Fab antibody fragments, and small proteins. It is useful for purifying albumin (K <sub>i</sub> =1.9 μM).
Targets(IC <sub>50</sub> )	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481