

DSPE-PEG3-azide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3-azide (DSPE-PEG3-N3) is a lipid-containing azide utilized in forming micelles for use as drug delivery nanoparticles. As a click chemistry reagent, it features an azide group capable of engaging in copper-catalyzed azide-alkyne cycloaddition (CuAAC) reactions with alkyne-containing molecules. Additionally, it can participate in strain-promoted azide-alkyne cycloaddition (SPAAC) reactions with molecules containing DBCO or BCN groups.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481