

6-Arm-PEG-FA (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	6-Arm-PEG-FA (MW 10000) [6-Arm-PEG-Folate (MW 10000)] is a PEG derivative modified with folic acid. Polyethylene glycol enhances solubility and stability, reduces immunogenicity of peptides and proteins, and inhibits nonspecific binding of charged molecules to modified surfaces. Folic acid (-FA) targets high-affinity folate receptors, facilitating receptor-mediated drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481