

## GMPCP trisodium

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |  |
|---------------|--|
| Description   | GMPCP (GpCp) trisodium (Compound 7 trisodium) is a guanine derivative and CD73 inhibitor, with $K_i$ values of 1110 nM for rat soluble CD73 and 410 nM for human soluble CD73. This compound is applicable in cancer research. |
| Targets(IC50) | CD73   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481