

4Arm-PEG-RB (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4Arm-PEG-RB (MW 2000) [4Arm-PEG-Rhodamine B (MW 2000)] is a four-armed fluorescent dye composed of Rhodamine B and PEG. It is suitable for applications in drug delivery and fluorescence imaging, with excitation and emission wavelengths of 546 nm and 610 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481