

BOC-NH-PEG-FA (MW 600)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BOC-NH-PEG-FA (MW 600) [BOC-NH-PEG-Folate (MW 600)] is a PEG derivative modified with Folic acid. The Folic acid (-FA) component has a high affinity for folate receptors, which makes it useful for targeting cell membrane receptors in drug delivery applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481