

RB-PEG-Mal (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	RB-PEG-Mal (MW 2000) (Rhodamine B-PEG-Mal (MW 2000)) is a fluorescent dye composed of Rhodamine B, PEG, and Maleimide (Mal). It is suitable for fluorescence labeling and imaging applications, with excitation and emission wavelengths of 546 nm and 610 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481