

Dabcyl-RGVVNASSRLA-Edans TFA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dabcyl-RGVVNASSRLA-Edans (TFA) [Dabcyl-CMV-Edans (TFA)] is a fluorescent substrate designed for human cytomegalovirus (HCMV) protease.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481