

## 4-Arm-PEG-OH (MW 5000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	4-Arm-PEG-OH (MW 5000) is a PEG derivative characterized by a four-arm structure and hydroxyl (OH) functional groups. It can react through the hydroxyl groups with functional groups on biomolecules such as proteins, peptides, and antibodies, making it suitable for drug delivery applications.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481