

NHS-PEG-CHO (MW 3400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NHS-PEG-CHO (MW 3400) (NHS-PEG-Aldehyde (MW 3400)) is a PEG derivative composed of an aldehyde group (-CHO), PEG units, and an NHS ester. The NHS ester can couple with amino acids or other amine-containing molecules. The aldehyde group is a reactive functional group that can interact with compounds containing amino or hydroxyl groups.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481