

## XIAP ligand-Linker Conjugate 1

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | XIAPligand-Linker Conjugate 1 is an E3 ubiquitin ligase ligand-linker conjugate that comprises the XIAPBIR2 ligand XB2M54 and a linker. XIAPligand-Linker Conjugate 1 is applicable in the synthesis of PROTAC ER $\alpha$ Degradar-1. |
| Targets(IC50) | E3 Ligase Ligand-Linker Conjugates,IAP   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481