

## Antibacterial agent 286

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Antibacterialagent 286 (compound 5) is an effective broad-spectrum antimicrobial agent that shows significant inhibitory activity against methicillin-resistant Staphylococcus aureus (MRSA) with a minimum inhibitory concentration (MIC) of 25 µg/mL.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481