

Antibacterial agent 286

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Antibacterial agent 286 (compound 5) is an effective broad-spectrum antimicrobial agent that shows significant inhibitory activity against methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) with a minimum inhibitory concentration (MIC) of 25 µg/mL. |
| Targets(IC50) | Antibacterial |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481