

RB-PEG-OH (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-OH (MW 10000) [Rhodamine B-PEG-OH (MW 10000)] is a fluorescent dye composed of Rhodamine B, PEG, and a hydroxyl group. It is useful for fluorescence labeling and imaging applications, with excitation and emission wavelengths of 546 nm and 610 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481