

4-Arm-PEG-OH (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Arm-PEG-OH (MW 2000) is a PEG derivative featuring a four-arm structure and hydroxyl (OH) functional groups. This compound can react with functional groups on biomolecules (such as proteins, polypeptides, antibodies, etc.) via the hydroxyl groups and is applicable in drug delivery systems.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481