

PDGFRA-IN-1

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PDGFRA-IN-1 (Compound 4p) is an inhibitor of platelet-derived growth factor receptor A (PDGFRA) with an IC50 of 1.25 $\mu$ M. It exhibits notable anti-cancer properties and is effective in eliminating cancer cells, including primary glioma cells derived from patients.
Targets(IC50)	PDGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481