

DSPE-PEG-Fmoc (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-Fmoc (MW 1000) is a PEG-derived lipid compound utilized for the formation of micelles, serving as nanoparticles for drug delivery. Fmoc acts as a widely used amino-protecting group, safeguarding amino groups from unintended reactions until it is removed for specific coupling reactions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481