

## Lipid OA2 dihydrochloride

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Lipid OA2 dihydrochloride is an ionizable cationic lipid used in the preparation of single-component lipid nanoparticles for siRNA delivery. This compound efficiently delivers SOCS1-siRNA to dendritic cells, silences the SOCS1 gene, and enhances antigen presentation and pro-inflammatory factor secretion abilities.
Targets(IC50)	Lipase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481